2024


MANSOURA COLLEGE

Mansoura High Institute Of Engineering And Technology

Chemistry Laboratory Manual


Practical Engineering Chemistry

For Engineering Chemistry


Common Laboratory Equipment


Contents

Experiment	page
(1) Titration of acids and bases	
Assignments (1)	
(2) Oxidation Reduction Titration	
Assignments (2)	
(3) Chloride Titration	
Assignments (3)	
(4) Determine the Equivalent weight of Iron by Chemical Displacement method.	
Assignments (4)	
(5) Estimation of CaO in Cement solution by rapid EDTA method	
Assignments (5)	
(6) Estimation of hardness of water by EDTA method	

COMMON LABORATORY GLASSWARES


Burette

Pipette

Test-tube


Graduated
cylinder

Measuring cylinder


Erlenmeyer Flask

Conical flask

Separating funnel

Volumetric flask


Beaker


Filter funnel


TITRATION ASSEMBLY

Correct method to note down the readings

Graduated Cylinder
The reading is 36.5 ml.

Burette
The reading is 27.8 ml

Constructions

A. GENERAL:

1. Keep your bags in the cupboards below the working table
2. First-aid kit is available for emergency use only in the laboratory. Band-aids for minor cuts are also available in the First-aid kit. Notify your instructor or the technicians if you use safety items.
3. Notify your instructor if any accidents and/or injuries, regardless of their severity. If you need medical treatment, you will be promptly taken to the nearby Health Center.
4. Learn the location and use of safety equipments, including the eyewash bottle, fire extinguisher, and sand bucket.
5. Work cautiously with chemicals only after you have learned about their potential hazards as well as the chemical properties. Laboratory has a catalogue of MSD (Material Safety Data Sheet) sheets that contain all the information about chemicals.
6. Wash your hands well before leaving the laboratory.
7. Keep your hands away from your face, while working.
8. Handle the apparatus and chemicals carefully.
9. Leave plenty of tap water after discarding the waste in the sink.
10. In the event of a chemical spill, large or small, consult your laboratory instructor or the technician as to the appropriate method of clean-up.

B. HANDLING OF CHEMICALS & WASTE DISPOSAL IN THE LABORATORY:

1. To avoid spattering of acids which can cause burns, always add acid to water. Never add water to acid.
2. Before taking any reagent, you must carefully read the label on the bottle. Many chemicals have similar names however they may exhibit different properties viz. concentration level, etc.
3. To avoid unnecessary waste, obtain only the required amount of chemicals in an experiment. Your instructor will tell you the proper procedure for dispensing liquids and solids.
4. Never return unused chemicals to the reagent bottle without prior permission of the instructor.
5. Follow scrupulously the instructor's instructions in case of disposing the chemicals. Dispose of non-hazardous, water soluble substances in the sink, and put insoluble materials such as filter paper in waste basket.
6. Broken glass must be put into the containers specified for that purpose.
7. Before leaving the Laboratory please ensure, clean off the surface. Remove matches & papers and wipe down the surface with wet cloth.

C. EYE PROTECTION:

1. If you get an irritating substance in your eye, move quickly to the eye washer, and wash your eyes thoroughly for at least 15 minutes. Do not take this incidence as a common one. Have someone notify the instructor of the accident so that you can be taken to the nearby Health Center immediately.

3. Remove contact lenses while performing experiment in the laboratory.

D. FIRE HAZARD:

1. In case of fire bring the fact immediately in the notice of concerned laboratory instructor.
2. Do not dry chemicals in a drying oven or heat any materials with an open flame unless specifically directed to do so by the laboratory instructor.

E. CONTACT & INGESTION HAZARD:

1. If you spill a corrosive substance on your skin or clothing, wash it off with plenty of water for 15 minutes. Notify the instructor of any spillage as soon as possible; he/she will provide any necessary secondary treatment and will arrange for your transportation to the Health Center, if necessary.
2. Never eat, drink, or taste anything in the laboratory.
3. Smoking & use of cell phones are strictly prohibited in the laboratory.

Experiment (1)

Titration of Acids and Bases

- **Objectives:**

- Standardize a sodium hydroxide solution
- Determine the molarity of an unknown hydrochloric acid solution
- Understand the use of indicators in titrations
-
- Learn proper pipetting technique
- Learn to titrate a strong acid with a strong base

- **Safety Notes:**

- Eye protection must be worn at all times.
- Hydrochloric acid and sodium hydroxide are caustic and should not come in contact with your skin or clothing.
- Wear gloves when handling these chemicals. A lab coat or lab apron is recommended.

- **Discussion:**

Titration describes a process where the concentration of an unknown substance is determined by comparing it with a solution of known concentration. The concept that makes titrations possible is finding the equivalence point, i.e., identifying when the quantity of the unknown substance is equal to the quantity of the known substance.

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The equivalence point is found in a titration by adding trace amounts of a substance, called an indicator, which turns color when the equivalence point is reached. When a strong acid is titrated with a strong base, or vice versa, the pH of the solution will be about 7.0 at the equivalence point. Phenolphthalein (**Ph.Ph**) is the indicator used in this experiment. Phenolphthalein (**Ph.Ph**) is colorless in acidic solutions and turns pink in alkaline solutions.

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This experiment will be done in two parts:

- (1) preparation and standardization of a 0.1 M sodium hydroxide solution and
- (2) determination of the acid concentration in an unknown sample.

Standardization is the process of determining concentration in an unknown solution by titrating it with a solution of known concentration.

Ionic substances dissociate completely when placed in water. Strong acids are ionic substances that form H^+ (or more correctly, H_3O^+) when placed in water and strong bases are ionic substances that form OH when dissolved in water.

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
Hydrochloric acid, HCl, is the strong acid used in this experiment and sodium hydroxide, NaOH, is the strong base. The approximate concentrations of our solutions will be 0.1 M, and the goal of the experiment is to determine the exact concentration of the unknown acid solution.

Acids contain hydrogen and dissociate in water to give H⁺:

Bases contain the hydroxyl group and dissociate in water to give OH⁻:


An acid (such as HCl) and a base (such as NaOH) react to form water and a salt. In this reaction, we say the acid and the base are **neutralized**:

Throughout this experiment the NaOH solution will be placed in the burette, and the hydrochloric acid solution will always be in the beaker. The phenolphthalein (**Ph.Ph**) indicator is colorless in acidic solution, and it will turn pink at the equivalence point. A small excess of OH⁻ is what causes the indicator to change from colorless to pink, which is the titration end point. (*Note: You should stop adding NaOH when the solution in the beaker holds a light pink color for about 10 seconds!*).

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• **Procedure:**

1. Rinse your burette twice with about 5 mL of the NaOH solution.

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2. Fill the burette with the NaOH solution and position it above an Erlenmeyer flask
(Note: be sure there are no air bubbles in the tip of the burette).

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
Figure 1. How to read a buret volume.

3. Obtain a 10.00 mL pipet and rinse it twice with about 5 mL of the unknown HCl acid solution.

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4. Transfer 10.00 mL of the unknown HCl acid solution into the Erlenmeyer flask.

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5. Add 1-2 drops of phenolphthalein indicator solution into the flask.

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6. Before you begin the titration, be sure to record the initial volume of NaOH in the burette.

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7. Begin the titration, swirling the solution in the flask as you add NaOH in a drop-wise approach (*Note: Your solution will initially turn pink, and then fade back to colorless when swirled. The pink color will remain longer as you approach the end point of the titration.*)

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8. Record the final volume of NaOH in the burette when the solution in the Erlenmeyer flask remains light pink color for about 10 seconds.

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9. Repeat this procedure two more times, refilling the NaOH in the burette as needed.

• **Results:**

Initial volume of NaOH (mL)	Final volume of NaOH (mL)	Total volume of NaOH used (mL)	Mean volume of NaOH used (mL)

Total volume NaOH used = Final volume NaOH – Initial volume NaOH

• **Calculations:**

Use the following relationship to calculate the HCl concentration:

$$V_a \times N_a = V_b \times N_b$$

where, V_a : The volume of NaOH solution, mL

N_a : The Normality of the NaOH solution

V_b : The volume of HCl solution, mL

N_b : The Normality of the HCl solution

This simplifies to give:

$$\frac{V_a N_a}{V_b} = N_b$$

Experiment (2)


Oxidation Reduction Titration

• Introduction:


Oxidation-reduction reactions (also known as redox reactions) are reactions that usually involve transfer of electrons. To determine the number of electrons transferred, oxidation states are assigned. Oxidation states of atoms are numbers that help chemists keep track of electrons during a reaction. Each atom in an equation can be assigned an oxidation state according to certain rules. If the oxidation state of an atom increases as you go from the reactants to the products in an equation, oxidation has occurred (electrons have been lost); if the oxidation state decreases, reduction has occurred (electrons have been gained).

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
In this experiment you will use a standard solution of potassium permanganate (KMnO_4) to determine the iron (as Fe^{2+}) in an unknown solution. Permanganate ion reduces to a manganese (II) ion in the acidic solution. This reaction requires 5 electrons and 8 hydrogen ions:

Only one electron is necessary to reduce Fe (III) to Fe (II)

Therefore, 1 mole of MnO_4^{-} (the oxidizing agent) reacts with 5 moles of Fe^{2+} (the reducing agent) to form 5 moles of Fe^{3+} and 1 mole of Mn^{2+} . Thus, in net ionic form:

The 1:5 mole ratio with respect to the amounts of MnO_4^{-} and Fe^{2+} consumed will provide the stoichiometric basis for all the calculations in this experiment.

The permanganate ion acts as its own indicator, as MnO_4^{-} is highly colored while Mn^{2+} is essentially colorless. The product of oxidation, the Fe^{3+} ion, is itself, slightly colored. To avoid any possible interference with the equivalence point determination a little phosphoric acid, H_3PO_4 , is added to complex Fe^{3+} to a completely, colorless ion.

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• **Procedure:**

1. Rinse the burette with about 3 mL of your standard KMnO_4 solution.

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2. Fill the burette using a plastic funnel. Clear all bubbles from the tip. Remove any drops from the tip.

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3. Record the initial volume, before starting each titration. The top of the meniscus will be read at the beginning and end of the titration.

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
4. Transfer 10.00 mL of the unknown Ferrous ammonium sulfate into the Erlenmeyer flask.

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5. Titrate with KMnO_4 until the appearance of a very faint pink color that persists for 30 seconds. *Note that as you approach this endpoint, the pink color will begin to persist for longer periods of time, before disappearing. At this point, it is advisable to add titrant slowly. You can even add partial drops from the burette tip - simply touch the flask to the droplet on the tip of the burette and then rinse the wall of the flask with a squirt of water from your wash bottle.*

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9. Repeat this procedure two more times, refilling the KMnO_4 in the burette as needed.

• **Results:**

Initial volume of KMnO ₄ (mL)	Final volume of KMnO ₄ (mL)	Total volume of KMnO ₄ used (mL)	Mean volume of KMnO ₄ used (mL)

Total volume KMnO₄ used = Final volume KMnO₄ – Initial volume KMnO₄

• **Calculations:**

Use the following relationship to calculate the Ferrous ammonium sulfate concentration:

$$V_a \times N_a = V_b \times N_b$$

where, V_a : The volume of KMnO₄ solution, mL

N_a : The Normality of the KMnO₄ solution

V_b : The volume of FAS solution, mL

N_b : The Normality of the FAS solution

This simplifies to give:


$$\frac{V_a N_a}{V_b} = N_b$$

Experiment (3)


Chloride Titration

• **Introduction:**

This method determines the chloride ion concentration of a solution by titration with silver nitrate. As the silver nitrate solution is slowly added, a precipitate of silver chloride forms.

The end point of the titration occurs when all the chloride ions are precipitated. Then additional silver ions react with the chromate ions of the indicator, potassium chromate, to form a red-brown precipitate of silver chromate.

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This method can be used to determine the chloride ion concentration of water samples from many sources such as seawater, stream water, river water and estuary water. Seawater is used as the example here.

• **Safety:**

- Lab coats, safety glasses and enclosed footwear must always be worn in the laboratory. The chromate solution needs to be prepared and used with care as

chromate is a known carcinogen. Silver nitrate solution causes staining of skin and fabric (chemical burns). Any spills should be rinsed with water immediately.

• **Equipment Needed:**

- burette and stand.
- 250 mL conical flasks.
- 10 mL measuring cylinders.

• **Procedure:**

1. Rinse the burette with about 3 mL of your standard AgNO_3 solution.

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3. Pipette a 10 mL aliquot of NaCl into a conical flask and add about 1 mL of chromate indicator.

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4. Titrate the sample with 0.1 mol L^{-1} silver nitrate solution. Although the silver chloride that forms is a white precipitate, the chromate indicator initially gives the cloudy solution a faint lemon-yellow color (**Figure 1**).

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Figure 1 Before the addition of any silver nitrate the chromate indicator gives the clear solution a lemon-yellow color.

The endpoint of the titration is identified as the first appearance of a red-brown color of silver chromate (**Figure 2**).

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Figure 2: **Left flask:** before the titration endpoint, addition of Ag^+ ions leads to formation of silver chloride precipitate, making the solution cloudy. The chromate indicator gives a faint lemon-yellow color. **Centre flask:** at the endpoint, all the Cl^- ions have precipitated. The slightest excess of Ag^+ precipitates with the chromate indicator giving a slight red brown coloration. **Right flask:** If addition of Ag^+ is continued past the endpoint, further silver chromate precipitate is formed and a stronger red brown color result. NB: The titration should be stopped

when the first trace of red brown color is observed. Using an incompletely titrated reference flask for comparison is a helpful way to identify the first appearance of red brown coloration.

5. Obtain the final volume reading from the calibration scale on the burette.

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9. Repeat this procedure two more times, refilling AgNO_3 in the burette as needed.

• **Results:**

Initial volume of AgNO_3 (mL)	Final volume of AgNO_3 (mL)	Total volume of AgNO_3 used (mL)	Mean volume of AgNO_3 used (mL)

Total volume AgNO_3 used = Final volume AgNO_3 – Initial volume AgNO_3

• **Calculations:**

Use the following relationship to calculate the Ferrous ammonium sulfate concentration:

$$V_a \times N_a = V_b \times N_b$$

where, V_a : The volume of AgNO_3 solution, mL

N_a : The Normality of the AgNO_3 solution

V_b : The volume of NaCl solution, mL

N_b : The Normality of the NaCl solution

This simplifies to give:

$$\frac{V_a N_a}{V_b} = N_b$$

Experiment (4)


Determine the Equivalent weight of Iron by Chemical Displacement method.

• **Object:**

To determine the Equivalent weight of Iron by Chemical Displacement method. The Equivalent weight of Copper is 63.5.




• **Theory:**

Iron displaces copper from a solution containing copper ions.

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
The quantity of copper sulphate remaining in solution, after the chemical displacement is estimated by Iodometric titration method.

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The iodine so liberated remains dissolved in excess of KI and is proportional to the amount of copper sulphate which is then titrated against sodium thiosulphate using starch as indicator. At the end point the blue color disappears and a white ppt. of cuprous iodide is obtained.

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Under chemical displacement conditions:

$$\frac{\text{Equivalent weight of copper}}{\text{Equivalent weight of Iron}(x)} = \frac{\text{Wt.of Cu.Deposited on Iron strip (C)}}{\text{Wt.of Iron goes into solution (D)}}$$

- **Chemicals:** N/10 $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$, N/10 sodium thiosulphate ($\text{Na}_2\text{S}_2\text{O}_3 \cdot 5\text{H}_2\text{O}$), solid KI, dilute H_2SO_4 .
- **Glassware:** Burette, Pipette, Beaker, conical flask, funnel, measuring cylinder.
- **Indicator:** Freshly prepared starch solution.
- **End point:** Disappearance of blue color.
- **Procedure:**

1. Take an iron strip (4 cm x 1 cm) and clean it with sandpaper. Weigh the cleaned iron strip accurately and place it in a clean 250 ml beaker. Pour 100 ml of CuSO_4 solution (N/10) into it and allow the strip to stand in beaker for about 30 minutes.

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2. With the help of forcep, carefully withdraw iron strip from the beaker and place it on a porcelain plate contained in a desiccator (using CaCl_2 as desiccators).

3. The quantity of copper sulphate remaining in solution after the chemical displacement is estimated by Iodometric titration method.

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4. The dried iron strip (containing the deposited copper) is then carefully weighed.

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5. Pipette out 25 ml of CuSO_4 solution (solution after chemical displacement). Now add 1 gm of KI in a conical flask, mix well and cover the mouth of conical flask with watch glass and allow the mixture to stand for 2 to 5 minutes. in the dark.

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6. Now titrate the liberated iodine with N/10 sodium thiosulphate solution. The brown color of iodine becomes fainter at the addition of sodium thiosulphate solution. When very light-yellow color remains add 5 drops of starch solution. It forms deep blue Iodine-starch complex Now add further Hypo solution drop by drop till blue color disappears. This is the end point.

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	Volume of CuSO ₄ solution taken (ml)	Initial	Final	sodium thiosulphate solution used (ml)

• **Observations:**

1. Initial weight of Iron strip, A = _____ gm
2. Weight of iron strip + copper after drying, B = _____ gm.
3. Weight of copper deposited on iron strip, C = gm.
4. Weight of iron which goes into solution (FeSO₄), D = A + C – B _____ gm.
5. Equivalent weight of copper = 63.5.

• **Calculations:**

1) Initial conc. of CuSO₄ solution = Normality x Equ.wt.
= 1 / 10 x 63.5 = 6.35 gm/lit.

2) $N_1 V_1 = N_2 V_2 \Rightarrow N_1 = \frac{N_2 V_2}{V_1}$

3) Final conc. of CuSO₄ = Normality x Equ.wt.

4) $\frac{\text{Equivalent weight of copper}}{\text{Equivalent weight of Iron(x)}} = \frac{\text{Wt.of Cu.Deposited on Iron strip (C)}}{\text{Wt.of Iron goes into solution (D)}}$

5) Percentage error =

• **Result:**

1. The Equivalent of Iron =
2. Percentage error =

Experiment (5)

Estimation of CaO in Cement Solution by rapid EDTA method

- **Aim:**

To determine the percentage of calcium oxide in the given sample of cement solution using Standard EDTA solution.

- **Principle:**

Cement contains compounds of calcium, aluminum, magnesium, iron, and insoluble silica. When dissolved in acid, silica remains undissolved. On treating with ammonia, aluminum and iron can be precipitated as their hydroxides and separated. The provided cement solution contains calcium and magnesium ions. The constituents of Portland cement are: - CaO (60-67%), SiO₂ (17-25%), Al₂O₃ (3-8%), Fe₂O₃ (0.5-6%), MgO(0.1 -4%), SO₃(1-3%), K₂O & Na₂O (0.5-1.5%) and CaSO₄(3-5%).

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To estimate the calcium content in the given solution, a known volume of cement solution in presence of Magnesium, Calcium ions is titrated with standard EDTA solution using Patton & Reader's indicator in the pH range 12-14. The indicator combines with Calcium ions to form a wine-red colored Calcium-Indicator

complex (Wine red). Near the end point, when free calcium ions are exhausted in the solution, further addition of EDTA, dissociates Calcium Indicator complex, consumes the calcium ions, and release free indicator which is blue in color. Therefore, color change is wine red to blue.

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• **Procedure:**

Part A: Preparation of standard solution of Disodium salt of EDTA

Weigh out the given EDTA crystals accurately into a 250 ml volumetric flask. Add quarter test tube of ammonia. Dissolve in distilled water and dilute up to the mark, mix well.

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Part B: Estimation of CaO in the given cement solution

Pipette out 25 cm³ of the given Cement solution into a clean conical flask. Add 5 ml of glycerol, 5 ml of diethyl amine and 10 ml of 4N NaOH solution. Add 3-4 drops of Patton and Reeder's indicator. Titrate against standard EDTA solution till the

color change from wine red to clear blue. Repeat the titration to get concordant values.

RESULT: CaO in the given sample of cement solution = -----

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• **Observation:**

Part A: Preparation of standard solution of Disodium salt of EDTA

Weight of weighing bottle + EDTA salt (W_1) = ----- g

Weight of empty weighing bottle (W_2) = ----- g

Weight of EDTA transferred ($W_1 - W_2$) = ----- g

$$\text{Molarity of EDTA} = \frac{\text{Weight of EDTA taken } (W_1 - W_2) \times 4}{\text{Molecular weight of EDTA (372)}} = \text{-----}$$

Part B: Estimation of CaO in the given cement solution

Burette reading	I	II	III
Final burette reading			
Initial burette reading			
Volume of EDTA run down (ml)			

• **Observation:**

◆ Volume of EDTA consumed by 25 cm³ of cement solution = -----

- ◆ Weight of Cement in 25ml =0.1g
- ◆ 1000 cm³ of 1 M EDTA = 56.08 g of CaO (Molecular mass of CaO =56.08)
- ◆ ----- cm³ of ----- M EDTA = $\frac{X \times 56.08}{1000 \times 1}$ g of CaO = -----
- ◆ 25 cm³ of the Cement solution contains ----- g of CaO
- ◆ % of CaO in the given sample of cement solution = $\frac{X \times 100}{0.1}$ =-----
- ◆ % of CaO in the given sample of cement solution = -----

Experiment (6)

Estimation of hardness of water by EDTA method

- **Aim:**

To estimate the total, permanent and temporary hardness of water by complexometric (EDTA) method.

- **APPARATUS:**


1. Conical flask
2. Burette
3. Pipette
4. Burette stand
5. Spatula

- **CHEMICALS:**

1. Buffers solution
2. Eriochrome black-T indicator (EBT)
3. EDTA Solution

- **PRINCIPLE:**


EDTA forms colorless, stable complexes with Ca^{2+} & Mg^{2+} ions present in water at pH 9-10. To maintain pH at 9-10, NH_4Cl & NH_4OH buffer is used. Eriochrome Black-T (EBT) is used as Indicator. Hard water sample with EBT indicator forms unstable, wine red colored complexes with Ca^{2+} & Mg^{2+} present in water.

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The metal indicator complex is less stable than metal EDTA complex. So, this metal indicator complex is titrated with std. EDTA solution. Now colorless metal-EDTA complex is formed by releasing EBT indicator which is blue color. So, the color change from wine red to blue indicates the end point.

(Colorless stable Indicator)(blue Complex in color)

One mole of EDTA reacts with one mole of complex, hence equimolar reaction. (n=1).

• **Procedure:**

Preparation of standard hard water:

Dissolve 1 gm of CaCO_3 (dry) in minimum (2-3 drops) quantity of dilute HCl and evaporate the solution on water bath. Dissolve this Solution in small amount of distilled water and then transfer into 1000 ml standard flask. Make up the solution up to the mark and shake well for uniform concentration.

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Standardization of EDTA Solution:

Pipette out 20 ml of std. hard water solution into a conical flask, add 2 ml of buffer solution and 2-3 drops of EBT indicator and titrate the wine red color complex solution with EDTA solution taken in the burette, till the color changes to blue color. Note the reading. Repeat the titration to get concurrent values.

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S.No.	Vol. of standard hard water	Burette reading		Vol. of EDTA required
		Initial	Final	

$$\frac{M_1 V_1}{n_1} = \frac{M_2 V_2}{n_2} \quad (\text{Note that } n_1 = n_2 = 1)$$

$$M_1 V_1 (\text{Hard water}) = M_2 V_2 (\text{EDTA solution})$$

• **Estimation of Total Hardness:**

Pipette out 50 ml of hard water sample (tap water) into 250 ml conical flask, add 2 ml of buffer solution and 3 drops of EBT-Indicator. Titrate the wine-red colored solution with EDTA, till the blue color is obtained. Note the reading. Repeat the titration to get concurrent values.

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S.No.	Vol. of standard hard water	Burette reading		Vol. of EDTA required
		Initial	Final	

$$\frac{M_3 V_3}{n_3} = \frac{M_2 V_2}{n_2} \quad (\text{Note that } n_3 = n_2 = 1)$$

$$M_3 V_3 (\text{Tap water}) = M_2 V_2 (\text{EDTA solution})$$

$$\text{Total Hardness} = M_3 \times 100 \times 1000 \text{ ppm}$$

• **Estimation of Permanent Hardness:**

Pipette out 100 ml of tap water and boil the solution till the volume is reduced to half i.e. 50 ml (All bicarbonates of Ca^{2+} & Mg^{2+} decomposes to CaCO_3 and $\text{Mg}(\text{OH})_2$ respectively). Cool the solution and filter the water into conical flask. Now

add 2 ml of buffer solution and 3 drops of EBT indicator, till the blue color is obtained. Note the readings. Repeat the titration to set concurrent values.

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S.No.	Vol. of standard hard water	Burette reading		Vol. of EDTA required
		Initial	Final	

$$\frac{M_4 V_4}{n_4} = \frac{M_2 V_2}{n_2} \quad (\text{Note that } n_4 = n_2 = 1)$$

$$M_4 V_4 \text{ (Heated tap water)} = M_2 V_2 \text{ (EDTA solution)}$$

$$\text{Permanent hardness} = M_4 \times 100 \times 1000 \text{ ppm}$$

$$\text{Temporary hardness} = \text{Total hardness} - \text{Permanent hardness}$$

• **RESULT:**

Temporary hardness of water = _____ ppm

Permanent hardness of water = _____ ppm

Total hardness of water = _____ ppm